

## **Audioscan introduces multiple software improvements – the most user-friendly Verifit ever!**

*Including dramatically faster data-entry and copy functions, enhanced fitting formulas and user sound file improvements.*

**Dorchester, Ontario – January 5, 2016** – Audioscan today introduced new user-friendly software for its market-leading Verifit, Axiom and RM500SL brands. The Verifit2 fully binaural model boasts more than nine new features including time-saving data enhancements where copying, deleting and entering audiometric data now need up to four times fewer clicks. All HL, BCT and UCL data are entered for both ears on just one single screen. Symmetric losses and RECD take half as long to enter with the new Audiogram/RECD copy function. In addition, DSL fitting has received a major upgrade with new sound field transducer types, target interpolation for incomplete audiograms, expanded 8kHz targets, up from 6Khz and a 10-point RECD, up from nine.

Audioscan president Jim Jonkman said, “This software release extends Audioscan’s famous ease-of use benefits in valuable ways for the hearing care professionals that we serve. I am proud of the work our world-class development team has delivered across all three product lines.”

Additional refinements in this new software release include:

1. tonal targets in NAL-NL2 (all models),
  2. monitor headphone level boost (Verifit2),
  3. selectable day-of-week for calibration expiry (all models),
  4. expanded dual table view in Speechmap (Verifit2),
  5. extended time for user sound files on USB (3x longer – up to 1 minute for all models),
  6. Korean speech passage (Verifit2, Axiom),
  7. new Shallow CIC instrument type (all models),
  8. ABR threshold entry for Axiom,
- and more.

Audioscan owners can download this powerful software for their specific model at no charge by clicking on <http://software.audioscan.com> or contact their local Audioscan sales representative for friendly assistance.

### **Press contact:**

John Pumford, Au.D.

Director of Audiology and Education

Audioscan

[john.pumford@audioscan.com](mailto:john.pumford@audioscan.com)

519-268-3313

For additional information, go to <http://audioscan.com>; <http://audiologyonline.com/expo/audioscan>

email: [info@audioscan.com](mailto:info@audioscan.com) Phone: 1-800-265-2093

©2017 Audioscan, a division of Etymonic Design Inc. All rights reserved.

Audioscan is the leading manufacturer of hearing instrument fitting systems in North America. The company manufactures the Verifit®2, Axiom®, and the RM500SL™ systems. Audioscan pioneered the Speechmap® sound visualization and the use of real speech for best-science hearing aid verification.